	STOOL COLLECTION FOR OCCULT BLOOD TEST	cod. S 51 date 11/04/2024 rev. 9
---	---	--

Dear Sir/Madam,

Please be advised that it is important to follow the instructions below to perform this test.

STOOL SAMPLE COLLECTION

The vials can be stored at room temperature.

- o Prepare a clean container to collect feces so that sampling can be done.
- o If 3-sample testing is required, it is necessary to collect each sample on 3 different stool releases.


Warnings

The preservative contains sodium azide.

- o Do not leave the container within reach of children.
 - o In case of skin contact, wash thoroughly with water.
 - o In case of contact with the eyes, wash thoroughly with water and contact the nearest emergency room.
 - o Do not ingest. In case of accidental ingestion, contact the nearest emergency room.
- In case of access to the emergency room, bring this preparation with you and show it to the medical staff.

STOOL SAMPLING

- o Collect the feces in a clean container. Be careful not to contaminate with urine or menstrual blood.
- o After collection, specimens should be stored in a refrigerator and delivered within 7 days during opening hours of the Collection Center.
- o **Caution:** do not spill the liquid from the vial.



Compiled by	Verified by	Approved by
Laboratory Director Dr. Laura Soldini	Quality Assurance Safety Environment Dr. Stefania Fazzone	Health Management Dr. Gabriele. Pellicciotta